

# Onix Fixed Ordering

## AspenLabs



**Prior to ordering, ensure all Clinical Records have been taken including ASAP Impressions, ASAP Scan, Scan Body Scan, and Headshots.**

**Step 1:** Navigate to the Patient Device Hub and find the "Fixed" device(s.) Click the checkbox with the proper device(s) and choose "Create Order."

**Step 2:** Click "Digital" then "Continue"

**Step 3:** Fill in the pre-liminary lab work as you normally would. *Remember, ASAP's lend the best outcome for the Onix Fixed Definitive cases.* Click "Save & continue"

**Step 4:** Ensure everything is correct as *you cannot go back.* Click "Looks good!"

**Step 5:** Proceed with filling out the script order form. *Note that "Lighthouse" is selected as the print location as Lighthouse will handle submitting the case to SKDLA on your behalf for design and manufacturing.*

**Step 6:** Monoblocks are recommended for all Onix Fixed cases. Please choose the arches that apply. Once you do that, a print location is required. If you are enabled to print in-office, you can choose "On-site," otherwise, Lighthouse/SKDLA will print & ship Monoblocks to your office.



# Onix Fixed Ordering

## AspenLabs



**Step 7:** Please ensure VDO, as well as tooth position and occlusal schema information is properly filled out. Once fully complete and necessary notes are added click "Save & Continue."

VDO\* 0 mm

Tooth position specifications\*

☐ Manually enter values ☒ Match existing denture

Midline 0 mm

Incisal edge 0 mm

High smile line 0 mm

Occlusal schema\*

☐ Lingualized (15 degree) ☐ Semi anatomic (22 degree)

☐ Combination (10 degree / 0 degree) ☒ Designer's preference

Cancel order NOTE: In order for this order to be moved to the design queue, it must be submitted first. To do so, once you have entered all the... Save draft & close Save & continue

UAT © 2025 - Aspen Dental

**Step 8:** Attach ASAP Scan, Scan Body Scan, and the 4 required Headshots (front, profile, smile, and high smile.) Then click "Continue."

AspenDental Clinical Account: 203312793 Patient: 0 - Hester AspenLabs - Patient device hub Chat for help Report Issues ID Schedule AspenLabs - 1-800-996-6470 UAT

Patient name: Hester Santiago Date of birth: 08/20/1958 Account number: 203312793 Digit: 203198015

Select the intraoral scans

New scan

July 10, 2025 [Recent] [Deselect set]

Upper [X] Side [X] Bite [X] Denture [X]

July 10, 2025 [Deselect set]

Upper [X] Side [X] Bite [X] Denture [X]

July 10, 2025 [Select set]

Upper [X] Side [X] Bite [X] Denture [X]

NOTE: In order for this order to be moved to the design queue, it must be submitted first. To do so, once you have entered all the required fields, please **continue** to the Order review screen and click "Submit".

Save draft & close Continue

UAT © 2025 - Aspen Dental (PM)

**Step 9:** Review the order and click "Submit"

AspenDental Clinical Account: 203312793 Patient: 0 - Lawrence AspenLabs - Patient device hub Chat for help Report Issues ID Schedule AspenLabs - 1-800-996-6470 UAT

Patient name: Lawrence Thornton Date of birth: 10/28/1959 Account number: 207807500 Digit: 207476188

Order #1464938 - Order review

Please verify all information is correct.

Upper arch [Signature]

Base shade: Original pink

Tooth shade: A1

Tooth size: Medium

Tooth shape: Soft

Lower arch [Signature]

Base shade: Original pink

Tooth shade: A1

Preliminary lab work

Existing denture: Both

ASAP impression: Both

Additional notes

Type notes here

Added: 10/20/2025 - 3:38 pm

Work instructions

VDO: No change

Midline: Match existing dentures

Incisal edge: Match existing dentures

High smile line: Match existing dentures

Occlusal schema: -

Tooth data

Missing: 1, 16, 17, 31, 32

Extractions: 18

Fixed restorations: None

Cancel order Submit

UAT © 2025 - Aspen Dental (PM)

For help scanning or ordering, please call the **Digital Help Desk at 1-800-996-6470** and follow the prompts for Digital Dentistry workflow training and Lighthouse Customer Service.

